



Moda[™]

Open up to brilliant sound

Four choices with performance and style.

Styled for Indigo[™] and Element[™]

unitron[®]
hearing 



Full performance. Comfortable fit. Brilliant sound.

Moda™, the stylish, fully featured 10A BTE is the most extensive offering of open fitting products available. Moda's brilliant sound quality is the result of a faster, more powerful digital signal processing (DSP) platform, along with premium innovations, adaptive and comfort features.

What makes Moda™ great?

- ▶ Available at four technology levels and four price points – from premium to entry
- ▶ Fully featured – Moda styled instruments boast the same set of high-value features found in the regular BTE styles of Indigo™, Element™ 16, Element™ 8 and Element™ 4
- ▶ Satisfy a wider range of hearing losses – choose from slim tube open fit or standard earmold





Now more clients can open up with Moda™

The open-fitting style of Moda is available in Indigo, Element 16, Element 8 and Element 4. That means no compromise in choice or performance. Moda styled instruments include more features and premium innovations than any other 10A BTEs in the premium, mid and entry level categories.

Indigo delivers control in a sophisticated, easy-to-fit, remarkable sounding hearing instrument. Indigo brings to the automatic world superior results in multiple listening situations.

Element 16 is an advanced, fully automatic hearing instrument delivering hands-free convenience combined with premium innovations and comfort features.

Element 8 provides hands-free performance featuring automatic directionality combined with premium innovations and comfort features.

Element 4 is an affordable hearing instrument featuring premium innovations and comfort features.

So comfortable and natural, it's easy to forget it's there.

A comfortable fit and natural listening experience help new wearers adjust to their hearing instruments and reduce listening fatigue, contributing to improved client acceptance and satisfaction.

- ▶ Slim tube open fittings do not occlude ears for a better, more natural sound and reduced listening fatigue
- ▶ Feedback cancellation improves listening comfort
- ▶ Very small 10A BTE fits comfortably behind the ear

Easy to fit. Easy to use.

- ▶ Same-day fitting and dispensing, with no impression or custom earmolds
- ▶ The same user-friendly U:fit™ fitting software is used with all Indigo and Element models
- ▶ As your clients' needs change, it's easy to switch between slim tube open fittings and earmolds
- ▶ Stylish open fitting 13 BTE, styled for Indigo, Element 16, Element 8 and Element 4 for longer battery life or direct audio input



indigo™
Moda

element™
Moda

Moda™ at-a-glance

	indigo™	e ¹⁶	e ⁸	e ⁴
Performance Profile				
Range of loss	Mild to severe	Mild to severe	Mild to severe	Mild to severe
Channels	16	16	8	4
Bands	16	16	8	8
Processing types	WDRC	WDRC Linear Limiting	WDRC Linear Limiting	WDRC Linear Limiting
Features				
Automatic programs	AutoPro4™	AutoPro2™	AutoMic™	
Manual listening programs	Up to 3	Up to 3	Up to 3	Up to 3
Ability to disable listening programs	• Manual programs only	•	•	•
Directional microphone system	16-band adaptive directional, directional, omni	Adaptive directional, directional, omni	Adaptive directional, directional, omni	Directional, omni
Data logging	•	•	•	•
AntiShock™		•	•	•
Speech enhancement		Speech enhancement LD	•	
Noise reduction	•	•	•	•
Wind noise manager	•	•	•	•
Feedback management	Adjustable phase canceller	Phase canceller	Phase canceller	Phase canceller
Automatic telephone solution	Easy-t	Easy-t	Easy-t	Easy-t
Slim Tube/Open fitting	•	•	•	•
Standard earmold	•	•	•	•
Program beep indicator	•	•	•	•
Low battery warning	•	•	•	•
Start up delay	•	•	•	•

Four products with full performance and style.