




## This is FUSE™

The first-ever crossover  
hearing instrument that puts  
the very best of open fit  
BTE inside the ear

# FUSE™

## Open evolves

Fuse™ is the shape-shifting new crossover hearing instrument that puts everything you love about open fit BTEs – a fully loaded, instant solution – out of sight and completely inside the ear. Fuse brings together all the most desired attributes and performance features of two hearing instrument categories into a completely new crossover category. Fuse continues the Unitron commitment to developing purpose-driven innovations that deliver true benefits to wearers and hearing healthcare professionals. Featuring the world’s first articulating joint, Fuse moves naturally within the ear canal for greater comfort and retention. And an entirely new take on venting delivers Fuse’s unique open, all inside fit.



OVERVIEW		Passport™ Fuse™	Latitude™ 16 Fuse™	Latitude™ 8 Fuse™
PERFORMANCE PROFILE	Max gain (2cc)	40 dB	40 dB	40 dB
	Max output (2cc)	109 dB	109 dB	109 dB
	Range of loss	mild to moderately-severe	mild to moderately-severe	mild to moderately-severe
	Battery size	10A	10A	10A
SIGNATURE FEATURES	SmartFocus™	•	•	
	Self learning	learning volume and smartFocus	learning volume and smartFocus	
	LearnNow™	•		
	Automatic program	autoPro™4	autoPro™3	autoPro™2
	Feedback management system	•	•	•
	Data logging	program, volume and smartFocus	program, volume and smartFocus	program and volume
ADDITIONAL FEATURES	Channels	20	16	8
	Manual listening programs	up to 3	up to 3	up to 3
	Microphone strategy	omnidirectional	omnidirectional	omnidirectional
	Speech enhancement LD	•	•	•
	Noise reduction	•	•	•
	AntiShock™	•	•	•
	Wind noise manager	•	•	•
	MyMusic™	•	•	•
ACCESSORIES	Smart Control (optional)	•	•	•



## Crossover to FUSE™

The modern shape of more satisfying listening

